

Arch Suspender Fitting Guide

You should following the Richie Brace® Standard Fitting Guide instructins first before applying the steps below for your Arch Suspender. *NOTE: The images below do not show the Richie Brace® inside of a shoe for demonstration purposes only and to better illustrate the use of the Arch Suspender. Your Richie Brace® should be fitted and adjusted when the brace is already inside your shoe.*

1. Richie Brace® is positioned and secured with front and back straps, according to Standard Brace instructions. Next, feed the Arch Suspender strap through the D-Ring on the upper limb support of the brace.
2. Pull upwards on the Arch Suspender strap while rolling your foot slightly through the outside. Feel the strap lift your arch.



3. When a comfortable lift support is felt under your arch, secure the strap back down upon itself across the top of the arch. Further tightening or loosening may be necessary after you begin walking with the device in place.

