

(Refer to www.richiebrace.com to assist with device selection)

QOL 1/10 Christine Place Capalaba QLD 4157 Ph 3823 1531 Fax 3823 1530 Email qoladmin@qol4feet.com.au	PODIATRIST & PATIENT INFORMATION		
	Practice Name _____		
	Podiatrist _____		
Address _____			
City: _____	State _____	Post Code _____	
Patient Name _____			
<input type="checkbox"/> Male <input type="checkbox"/> Female Age: _____ Height: _____ Weight: _____ Shoe Size: _____ Shoe Type _____			
Cast enclosed for <input type="checkbox"/> Left <input type="checkbox"/> Right			
PLEASE MARK MEDIAL AND LATERAL MALLEOLI ON NEGATIVE CAST!			

CLINICAL INFORMATION	
DIAGNOSIS:	Stance Evaluation
	Calcaneus alignment to leg: _____° inverted or _____° everted
	Leg alignment to floor: _____° varum or _____° valgum

RICHIE BRACE® PRESCRIPTION	
<input type="checkbox"/> RICHIE BRACE® (standard): Full Flexion Ankle Hinge Pivot. Can include enhancements for Posterior Tibial Tendon Dysfunction (select any or all): Medial Heel Skive <input type="checkbox"/> 2mm <input type="checkbox"/> 4mm <input type="checkbox"/> 6mm Navicular Accommodation <input type="checkbox"/> (please mark negative cast) Adjust Limb Uprights for Tibial Varum <input type="checkbox"/> Yes <input type="checkbox"/> No (see measurements above) FOR SEVERE PTTD, RECOMMENDED MEDIAL ARCH SUSPENDER (SEE BELOW)	
<input type="checkbox"/> RICHIE BRACE® RESTRICTED ANKLE PIVOT: Limits ankle motion, yet allows smooth contact phase of gait. Indications: DJD ankle & STJ, tarsal coalition, mild Charcot, lateral ankle instability, peroneal tendinopathy. ENHANCEMENTS (optional): <input type="checkbox"/> MEDIAL ARCH SUSPENDER – Adjustable lifting strap under talo-navicular joint for severe PTTD <input type="checkbox"/> LATERAL ARCH SUSPENDER – Adjustable lifting strap under calcaneal-cuboid joint for peroneal tendinopathy and severe lateral ankle instability.	
<input type="checkbox"/> RICHIE BRACE® DYNAMIC ASSIST: Full flexion pivot with spring hinges for dorsiflexion assist. Patient requirements: 1. Dropfoot 2. Ankle dorsiflexion to at least 90° to leg 3. Stable knee (must have all 3) Dynamic Assist covers are extended to the sulcus and the extension is reinforced to increase stiffness	
ALL RICHIE BRACES® HAVE THE FOLLOWING STANDARD FEATURES:	
<ul style="list-style-type: none"> ◆ Top Cover – EVA 190 3mm Blue/Green ◆ Brace Color –Black ◆ Heel Cup – Deep 	<ul style="list-style-type: none"> ◆ Cover Length – Met heads ◆ Orthotic Foot Plate – Intrinsic Balance to Perpendicular
<ul style="list-style-type: none"> ◆ Limb Uprights Supports – Aligned Perpendicular to Foot Plate ◆ Heel Stabilizer Bar - Included 	

Your prescription is now complete, unless you wish to make any modifications:

RICHIE BRACE® MODIFICATIONS			
NOTE: NON-STANDARD BRACE MODIFICATIONS MAY HAVE EXTRA CHARGES – SEE PRICING SHEET			
Top Cover	Length	Heel Cup	Medial Heel Skive
<input type="checkbox"/> Multiform 3mm red <input type="checkbox"/> black <input type="checkbox"/> blue <input type="checkbox"/> multicolour <input type="checkbox"/> <input type="checkbox"/> Diabetic (Plastazote/Poron)	<input type="checkbox"/> Standard <input type="checkbox"/> to Sulcus <input type="checkbox"/> to Toes <input type="checkbox"/> add poron cushion to extension	<input type="checkbox"/> 10 mm <input type="checkbox"/> 14 mm <input type="checkbox"/> 18 mm <input type="checkbox"/> Standard	For severe pronation control <input type="checkbox"/> 2mm <input type="checkbox"/> 4mm <input type="checkbox"/> 6mm

CAST AND ORTHOTIC MODIFICATIONS		
<input type="checkbox"/> Add Medial Arch Flange <input type="checkbox"/> Add Lateral Flange <input type="checkbox"/> other _____	<input type="checkbox"/> Orthotic shell Accommodation (please mark on cast and illustration below) <input type="checkbox"/> Navicular <input type="checkbox"/> Plantar Fascia _____mm <input type="checkbox"/> Styloid 5 th Met <input type="checkbox"/> Other:	<input type="checkbox"/> ¼ length heel raise _____mm (EVA 350 heel to sulcus tapered wedge) Supplied and billed as a separate item
SPECIAL INSTRUCTIONS:		Accommodation location(s): (mark on illustration and on cast)
		Cast Return <input type="checkbox"/>